THUASNE

VENOFLEX *Incognito* **ABSOLU**

VENOFLEX *Kokoon* **ABSOLU**

VENOFLEX *Simply* **COTON FIN**

VENOFLEX *Fantaisie*

VENOFLEX *Secret*

VENOFLEX *Secret* **OPAQUE**

VENOFLEX *City Confort*

VENOFLEX *City Confort* **COTON**

VENOFLEX *City Confort* **FIL D'ECOSSE**

VENOFLEX *Élégance*

VENOFLEX *Fast* **LIN**

VENOFLEX *Fast* **COTON**

VENOFLEX *Fast* **LAINE**

VENOFLEX *Fast* **AIR**

VENOFLEX *Soft & Care*

VENOFLEX *Clinic*

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Conserver cette notice.

en

MEDICAL COMPRESSION SOCKS, STOCKINGS OR TIGHTS FOR WOMEN AND MEN

Description/Destination

VENOFLEX devices are intended for use on the lower limbs, by adults with the indications presented hereafter and whose measurements correspond to the size chart.

Composition: polyamide - elastane
Components specific to certain product ranges (see label): viscose - linen - wool - cotton - silicone*

* Silicone for the stockings

Properties/Mode of action

VENOFLEX devices provide compression by applying pressure to the limb, thereby promoting venous return.

Indications

- 10-15 mmHg/Class 1: heavy legs, restless legs, evening oedema, telangiectasia or reticular veins (< 3 mm), prevention of functional disorders linked to the prolonged standing, prevention of venous thrombosis during long journeys, pregnancy and post-partum (prevention and treatment of minor chronic venous insufficiency and deep venous thrombosis prevention), orthostatic hypotension.
- 15-20 mmHg/Class 2: heavy legs, restless legs, evening oedema, telangiectasia or reticular veins (< 3 mm), following sclerotherapy for telangiectasia or reticular veins (< 3 mm), following sclerotherapy or surgery for varicose veins (≥ 3 mm), varicose veins (≥ 3 mm), prevention of venous thrombosis in pregnant women or post-partum, prevention of venous

thrombosis during long journeys, treatment of lymphoedema (maintenance phase).

- 20-36 mmHg/Class 3: varicose veins (≥ 3 mm), following sclerotherapy or surgery for varicose veins (≥ 3 mm), treatment of chronic venous oedema, trophic disorders (pigmentation, eczema, lipodermatosclerosis or chronic atrophie blanche), healed venous ulcer, chronic venous disease associated with pregnancy or the post-partum period, prevention and/or treatment of venous thrombosis (lower limb) and post-thrombotic syndrome, treatment of lymphoedema (maintenance phase).

Contraindications

Do not use in case of major dermatological conditions on an affected limb. Do not use in the event of known allergy to any of the components. Do not use in case of lower limb peripheral arterial disease (PAD) with ankle-brachial index (ABI) < 0.6 . Do not use in the event of decompensated heart failure. Do not use in the event of advanced diabetic microangiopathy (for a compression > 30 mmHg). Do not use in case of severe peripheral neuropathy of the affected limb. Do not use in the event of plegmasia cerulea dolens (painful blue phlebitis with arterial compression) of an affected limb. Do not use in case of extra-anatomical bypass of the affected limb. Do not use in the event of septic thrombosis.

Precautions

Strictly comply with your healthcare professional's prescription and recommendations for use. The risk/benefit ratio and the adequate level of pressure must be regularly re-assessed by a healthcare professional in the event of:

- PAD with ABPI between 0.6 and 0.9
- Advanced peripheral neuropathy
- Weeping or eczematous dermatitis.

In the event of discomfort, significant restriction, pain, variation in limb volume, deterioration of skin condition, infection, unusual sensations, change in the colour of the extremities, or changes in the performance of the product, remove the

product and seek the advice of a healthcare professional. Do not apply the product directly to an open wound without a dressing. Check that the anti-slip system can be applied to healthy, unbroken skin. Take care not to damage the device when fitting it, particularly with your nails. For hygiene, security and performance reasons, do not re-use the product for another patient. These devices are designed for a daytime use. Do not use when swimming. Do not use products on the skin (creams, ointments, oils, gels, patches...) before applying the device, it could damage the product. For some pathologies (or situations) such as venous thrombosis, the device must be used in combination with an anticoagulant treatment; consult a healthcare professional for advice.

Venoflex Clinic: These products do not contain any metal wires liable to interfere with a magnetic field. In the event of use of an electro-surgical knife, it is recommended that a sock be worn rather than a stocking in order to keep the device away from the knife. The sock must be far enough away to ensure it does not come into contact with the plate or electrodes of the electro-surgical knife.

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity or cutaneous dryness. Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

Choose the appropriate size to fit the patient, referring to the size chart.

Verify the product's integrity before every use. It is recommended that the device be washed before it is used for the first time. Handle the product with care when putting on (pay attention to long nails and jewellery).

How to put on VENOFLEX socks

① (A) Place your hand inside the sock and pinch the heel

① (B) Slip your foot into the sock until the heel is correctly positioned

① (C) Pull up to below the knee and do not turn the sock on itself

① (D) Rub your leg to spread the sock evenly
How to put on VENOFLEX FAST and SOFT & CARE socks?

② (A) Turn back the top of the sock to the heel

② (B) Slip your foot into the sock until the toes and heel are correctly positioned

② (C) Unroll the sock from the ankle to below the knee and do not turn the sock on itself

② (D) Rub your leg to spread the sock evenly
How to put on VENOFLEX stockings or tights

③ (A) Place your hand inside the stocking or tights and pinch the heel

③ (B) Slip your foot into the stocking or tights until the heel is correctly positioned

③ (C) Carefully pull up to the top of the thigh (do not pull too hard on the anti-slip band), or to the waist for tights

③ (D) Rub your leg to spread the stocking or tights evenly

Obtaining the adequate pressure may require overlapping/combining with other compression devices.

Care: machine washable at 40°C (normal or delicate programme for Incognito absolu and Fast laine). If possible use a washing net. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

Storage: store at room temperature, preferably in the original packaging.

Disposal: Dispose of in accordance with local regulations.

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©Thuasne - 2043101 (2022-11)
THUASNE UK Ltd